

# Recombinant Human Serum Amyloid A protein ab87757

## 1 图像

### 描述

产品名称	重组人Serum Amyloid A蛋白
纯度	> 95 % SDS-PAGE. ab87757 is purified using conventional chromatography techniques.
表达系统	Escherichia coli
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	<b>MGSSHHHHH SSGLVPRGSH MRSFFSFLGE</b> <b>AFDGARDMWR AYSDMREANY IGSDKYFHAR</b> <b>GNYDAAKRGV GGVWAAEAI S DARENIQRFF</b> <b>GHGAEDSLAD QAANEWGRSG KDPNHFRPAG LPEKY</b>

### 技术指标

Our **Abpromise guarantee** covers the use of **ab87757** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用 SDS-PAGE

形式 Liquid

### 制备和贮存

稳定性和存储 Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

pH: 8.00


Constituents: 0.242% Tris, 10% Glycerol (glycerin, glycerine)

### 常规信息

相关性 Function: Major acute phase reactant. Apolipoprotein of the HDL complex. Tissue specificity: Expressed by the liver; secreted in plasma. Disease: Note=Reactive, secondary amyloidosis is

characterized by the extracellular accumulation in various tissues of the SAA protein. These deposits are highly insoluble and resistant to proteolysis; they disrupt tissue structure and compromise function. Note=Elevated serum SAA protein levels may be associated with lung cancer. Similarity: Belongs to the SAA family. PTM: This protein is the precursor of amyloid protein A, which is formed by the removal of approximately 24 residues from the C-terminal end.

## 图片



15% SDS-PAGE analysis of 3µg ab87757.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors