

Recombinant Human SMP30 protein (denatured) ab113154

1 图像

描述

产品名称	重组人SMP30蛋白(denatured)
纯度	> 85 % SDS-PAGE.
表达系统	Escherichia coli
Accession	<u>Q15493</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
预测分子量	35 kDa including tags
氨基酸	1 to 299
标签	His tag N-Terminus

技术指标

Our **Abpromise guarantee** covers the use of **ab113154** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE
形式	Liquid
补充说明	This product was previously labelled as Regucalcin.


制备和贮存

稳定性和存储	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. pH: 8.00 Constituents: 1.2% Urea, 0.32% Tris HCl, 20% Glycerol (glycerin, glycerine)
--------	---

常规信息

功能	Gluconolactonase with low activity towards other sugar lactones, including gulonolactone and galactonolactone. Can also hydrolyze diisopropyl phosphorofluoridate and phenylacetate (in vitro). Calcium-binding protein. Modulates Ca(2+) signaling, and Ca(2+)-dependent cellular processes and enzyme activities.
序列相似性	Belongs to the SMP-30/CGR1 family.
细胞定位	Cytoplasm.

图片



15% SDS-PAGE showing ab113154 (3 µg) at approximately 35.4 kDa.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors