

Recombinant Human Stanniocalcin 1/STC protein ab63280

1 图像

描述

产品名称	重组人Stanniocalcin 1/STC蛋白
纯度	> 95 % SDS-PAGE. Please filter the product by an appropriate sterile filter before using it in the cell culture.
表达系统	HEK 293 cells
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	<pre> THEAEQNDSV SPRKSRVAAQ NSAEVVRCLN SALQVGCRAF ACLENSTCDT DGMVDICKSF LYSAAKFDTQ GKAFVKESLK CIANGVTSKV FLAIRRCSTF QRMIAEVQEE CYSKLNCSI AKRNPEAITE VVQLPNHFSN RYYNRLVRSI LECEDETVST IRDSLMEKIG PNMASLFHIL QTDHCAQTHP RADFNRRRTN EPQKLKVLRL NLRGEEDSPS HIKRTSHESA - ASDYKDDDDK </pre>

技术指标

Our **Abpromise guarantee** covers the use of **ab63280** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Western blot
	ELISA
	SDS-PAGE
形式	Lyophilized
补充说明	Previously labelled as Stanniocalcin 1.

制备和贮存

稳定性和存储	Shipped at 4°C. Upon delivery aliquot. Store at -80°C. Avoid freeze / thaw cycle. Constituents: 0.242% Tris, 0.116% Sodium chloride
--------	--

复溶

Add deionized water to prepare a working stock solution of approximately 0.5 mg/mL and let the lyophilized pellet dissolve completely.

常规信息


功能 Stimulates renal phosphate reabsorption, and could therefore prevent hypercalcemia.

组织特异性 Expressed in most tissues, with the highest levels in ovary, prostate, heart, kidney and thyroid. In the kidney, expression is confined to the nephron, specifically in the distal convoluted tubule and in the collecting tubule. Not detected in the brain, liver, spleen, peripheral blood leukocytes and adrenal medulla.

序列相似性 Belongs to the stanniocalcin family.

细胞定位 Secreted.

图片



LANE 1. M.W. marker - 14, 21, 31, 45, 66, 97 kDa. LANE 2. reduced and heated sample, 5µg/lane. LANE 3. non-reduced and non-heated sample, 5µg/lane

SDS-PAGE - Recombinant Human Stanniocalcin
1/STC protein (ab63280)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors