

Recombinant Human Stanniocalcin 2/STC-2 protein ab63281

1 图像

描述	
产品名称	重组人Stanniocalcin 2/STC-2蛋白
纯度	> 90 % SDS-PAGE. Please filter the product by an appropriate sterile filter before using it in cell culture.
内毒素水平	< 0.100 Eu/μg
表达系统	HEK 293 cells
Accession	<u>Q76061</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	TDATNPPEGPQDRSSQQKGRLSLQNTAEIQHCLVNAGDVGCG VFECFENN SCEIRGLHGICMTFLHNAGKFDAQGKSFIKDALKCKAHALRH RFGCISRK CPAIREMVSQLQRECYLKHDLCAAAQENTRVIVEMIHFKDLL LHEPYVDL VNLLTTCGEEVKEAITHSVQVQCEQNWGSLCSILSFCTSAIQ KPPTAPPE RQPQVDRTKLSRAHHGEAGHHLPEPSSRETGRGAKGERGSKS HPNAHARG RVGGLGAQGPSGSSEWEDEQSEYSDIRRAAADYKDDDDK
预测分子量	33 kDa
氨基酸	25 to 302
额外的序列信息	C-Terminal Flag-tag

技术指标	
Our Abpromise guarantee covers the use of ab63281 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	Western blot ELISA

形式	Lyophilized
补充说明	Previously labelled as Stanniocalcin 2.
制备和贮存	
稳定性和存储	Shipped at 4°C. Upon delivery aliquot. Store at -80°C. Avoid freeze / thaw cycle. Constituents: 0.242% Tris, 0.116% Sodium chloride
复溶	Add deionized water to prepare a working stock solution of approximately 0.5 mg/mL and let the lyophilized pellet dissolve completely.
常规信息	
功能	Has an anti-hypocalcemic action on calcium and phosphate homeostasis.
组织特异性	Expressed in a variety of tissues including muscle, heart, pancreas, kidney, spleen, prostate, small intestine, colon and peripheral blood leukocytes.
序列相似性	Belongs to the stanniocalcin family.
细胞定位	Secreted.

图片



LANE 1. M.W. marker - 14, 21, 31, 45, 66, 97 kDa. LANE 2. reduced and heated sample, 5µg/lane. LANE 3. non-reduced and non-heated sample, 5µg/lane.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors