

Recombinant Human TCR alpha protein ab112358

1 图像

描述	
产品名称	重组人TCR alpha蛋白
表达系统	Wheat germ
Accession	<u>Q6IPR8</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MAMLLGASVLILWLQPDWVNSQQKNDDQVKQNSPSLSVQEG RISILNCD YTNSMFDYFLWYKKYPAEGPTFLISSIKDKNEDGRFTVFL NKSAKHLS LHIVPSQPGDSAVYFCAVRYNNNDMRFGAGTRLTVKPNIQNP DPAVYQLR DSKSSDKSVCLFTDFDSQTNVSQSKSDSVYITDKTVLDMRSM DFKSNSAV AWSNKSDFACANAFNNSIIPEDTFFPSPESSCDVKLVEKSFE TDTNLNFQ NLSVIGFRILLKLVAGFNLLMTLRLWSS
预测分子量	57 kDa including tags
氨基酸	1 to 278
标签	GST tag N-Terminus
技术指标	
Our Abpromise guarantee covers the use of ab112358 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	SDS-PAGE Western blot ELISA
形式	Liquid
制备和贮存	

稳定性和存储

Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.31% Glutathione, 0.79% Tris HCl

常规信息


相关性

TCR (T cell antigen receptor) recognizes foreign antigens and translates such recognition events into intracellular signals that elicit a change in the cell from a dormant to an activated state. TCR is a heterodimer composed of either alpha and beta or gamma and delta chains. CD3 chains and the CD4 or CD8 coreceptors are also required for efficient signal transduction through the TCR.

细胞定位

Plasma membrane

图片



ab112358 on a 12.5% SDS-PAGE Stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors