

Recombinant human TGF beta 3 protein (Active) ab52313

[1 References](#) [2 图像](#)

描述


产品名称	重组人TGF beta 3蛋白(Active)
生物活性	ab52313 is fully biologically active. The ED ₅₀ was determined by cytotoxicity assay using the WHO standard 98/608 as a direct comparison is < 0.05 ng/ml, corresponding to a specific activity of > 2 x 10 ⁸ units/mg.
纯度	> 98 % SDS-PAGE. >98% HPLC analyses
表达系统	Escherichia coli
Accession	<u>P10600</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	ALD ^T NYCFRN LEENCCVRPL YIDFRQDLGW K ^W VHEPKGY ^Y ANFCSGPCPY LRSADTTHTS ^T VLGLYNTLNP EASASPCCVP QDLEPLTILY YVGRTPKVEQ LSMVVK ^S CK CS
预测分子量	26 kDa
氨基酸	301 to 412
额外的序列信息	(Gene symbol: TGFB3) (Gene ID: 7043)

技术指标

Our **Abpromise guarantee** covers the use of **ab52313** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	HPLC Western blot Functional Studies SDS-PAGE
形式	Liquid
补充说明	This product is manufactured by BioVision, an Abcam company and was previously called 4344



Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors