

# Recombinant Human TSC-1 protein ab87891

## 1 图像

### 描述

---

|       |  |
|-------|--|
| 产品名称  | 重组人TSC-1蛋白   |
| 纯度    | > 90 % SDS-PAGE.<br>ab87891 was purified by conventional chromatography, after refolding of the isolated inclusion bodies in a renaturation buffer. Endotoxin level: < 1.0 EU per 1 microgram of protein (determined by LAL method). |
| 表达系统  | Escherichia coli   |
| 蛋白长度  | Full length protein  |
| 无动物成分 | No   |
| 性质    | Recombinant  |
| 种属    | Human  |

### 技术指标

---

Our **Abpromise guarantee** covers the use of **ab87891** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

|      |                                   |
|------|-----------------------------------|
| 应用   | SDS-PAGE                          |
| 形式   | Liquid                            |
| 补充说明 | Previously labelled as Eotaxin 3. |

### 制备和贮存

---

|        |  |
|--------|--|
| 稳定性和存储 | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.<br>pH: 7.40<br>Constituents: PBS, 10% Glycerol (glycerin, glycerine) |
|--------|--|

### 常规信息

---

|       |  |
|-------|--|
| 功能    | Chemotactic for eosinophils and basophils. Binds to CCR3.                          |
| 组织特异性 | Ubiquitously expressed at low levels in various tissues including heart and ovary. |


**序列相似性**

Belongs to the intercrine beta (chemokine CC) family.

**细胞定位**

Secreted.

---

**图片**

15% SDS-PAGE showing ab87891 (3 $\mu$ g).

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

**Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

**Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors