

Recombinant Human Ube2L3/UBCH7 protein ab48775

1 图像

描述	
产品名称	重组人Ube2L3/UBCH7蛋白
纯度	> 95 % SDS-PAGE. Ube2L3/UBCH7 was overexpressed in E.coli and purified by using conventional chromatography techniques.
表达系统	Escherichia coli
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MAASRRLMKE LEEIRKCGMK NFRNIQVDEA NLLTWQGLIV PDNPPYDKGA FRIEINFPAE YPFKPPKITF KTKIYHPNID EKGQVCLPVI SAENWKPATK TDQVIQSLIA LVNDPQPEHP LRADLAE EYS KDRKKFCKNA EEFTKKYGEK RPVD
氨基酸	1 to 154

技术指标	
Our Abpromise guarantee covers the use of ab48775 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	SDS-PAGE
形式	Liquid
补充说明	This product was previously labelled as Ube2L3

制备和贮存	
稳定性和存储	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. pH: 7.4 Constituents: 0.0154% DTT, 1.19% HEPES, 10% Glycerol (glycerin, glycerine), 0.87% Sodium chloride

常规信息

功能	Accepts ubiquitin from the E1 complex and catalyzes its covalent attachment to other proteins. In vitro catalyzes 'Lys-11'-linked polyubiquitination. Involved in the selective degradation of short-lived and abnormal proteins. Down-regulated during the S-phase it is involved in progression through the cell cycle. Regulates nuclear hormone receptors transcriptional activity. May play a role in myelopoiesis.
组织特异性	Ubiquitous, with highest expression in testis.
通路	Protein modification; protein ubiquitination.
序列相似性	Belongs to the ubiquitin-conjugating enzyme family.
翻译后修饰	Ubiquitinated. The alteration of UBE2L3 protein levels during the S-phase of the cell cycle is due to ubiquitin-dependent proteasomal degradation.
细胞定位	Nucleus. Cytoplasm.

图片



14% SDS-PAGE gel loaded with ab48775 recombinant human Ube2L3/UBCH7 protein predicted molecular weight 17.9KDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

-
- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors