

# Recombinant human ULK2 protein ab125622

## 5 图像

### 描述

---

产品名称	重组人ULK2蛋白
生物活性	The specific activity of ab125622 was determined to be 102 nmol/min/mg.
纯度	> 80 % Densitometry. Affinity purified.
表达系统	Baculovirus infected Sf9 cells
Accession	<a href="#">Q8IYT8</a>
蛋白长度	Protein fragment
无动物成分	No
性质	Recombinant
种属	Human
预测分子量	98 kDa including tags
氨基酸	1 to 478

### 技术指标

---

Our [Abpromise guarantee](#) covers the use of **ab125622** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Western blot Functional Studies SDS-PAGE
形式	Liquid
补充说明	<a href="#">ab64311</a> (Myelin Basic Protein protein) can be utilized as a substrate for assessing kinase activity

### 制备和贮存

---

稳定性和存储	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 7.50 Constituents: 0.004% EGTA, 0.002% PMSF, 0.004% DTT, 0.79% Tris HCl, 0.003% EDTA, 25% Glycerol (glycerin, glycerine), 0.88% Sodium chloride
--------	---

This product is an active protein and may elicit a biological response in vivo, handle with caution.

## 常规信息

### 序列相似性


Belongs to the protein kinase superfamily. Ser/Thr protein kinase family. APG1/unc-51/ULK1 subfamily.

Contains 1 protein kinase domain.

### 翻译后修饰


Autophosphorylated.

## 图片




Functional Studies - Recombinant human ULK2 protein (ab125622)

The specific activity of ULK2 (ab125622) was determined to be 110.4 nmol/min/mg as per activity assay protocol




SDS-PAGE - Recombinant human ULK2 protein (ab125622)

SDS PAGE analysis of ab125622




SDS-PAGE - Recombinant human ULK2 protein (ab125622)

SDS PAGE analysis of ab125622



Functional Studies - Recombinant human ULK2 protein (ab125622)

Sample Kinase Assay showing the specific activity of ab125622 as 102 nmol/min/mg.



SDS-PAGE - Recombinant human ULK2 protein (ab125622)

SDS-PAGE analysis of ab125622.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors