

# Recombinant Human ULK4 protein ab125657

## 1 图像

### 描述

---

产品名称	重组人ULK4蛋白
纯度	> 75 % Densitometry. Purity was determined to be >75% by densitometry. Affinity purified.
表达系统	Baculovirus infected Sf9 cells
Accession	<b><u>Q96C45</u></b>
蛋白长度	Protein fragment
无动物成分	No
性质	Recombinant
种属	Human
预测分子量	97 kDa including tags
氨基酸	1 to 616

### 技术指标

---

Our **Abpromise guarantee** covers the use of **ab125657** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Western blot SDS-PAGE
形式	Liquid

### 制备和贮存

---

稳定性和存储	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 7.50 Constituents: 0.31% Glutathione, 0.002% PMSF, 0.004% DTT, 0.79% Tris HCl, 0.003% EDTA, 25% Glycerol (glycerin, glycerine), 0.88% Sodium chloride
--------	---

### 常规信息

---


序列相似性	Belongs to the protein kinase superfamily. Ser/Thr protein kinase family. APG1/unc-51/ULK1
-------	--

subfamily.  
Contains 5 HEAT repeats.  
Contains 1 protein kinase domain.

## 结构域

The protein kinase domain is predicted to be catalytically inactive.

## 图片



SDS-PAGE analysis of ab125657.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors