

Recombinant Human VILIP1 protein ab63276

1 图像

描述	
产品名称	重组人VILIP1蛋白
纯度	> 95 % SDS-PAGE.
表达系统	Escherichia coli
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	<div><div>MKHHHHHAS</div><div>MGKQNSKLAP</div><div>EVMEDLVKST</div><div>EFNEHELKQW</div><div>YKGFLKDCPS</div><div>GRLNLEEFQQ</div><div>LYVKFFPYGD</div><div>ASKFAQHAFR</div><div>TFDKNGDGTI</div><div>DFREFICALS</div><div>ITSRGSFEQK</div><div>LNWAFNMYDL</div><div>DGDGKITRVE</div><div>MLEIIEAIYK</div><div>MVGTVIMMKM</div><div>NEDGLTPEQR</div><div>VDKIFSKMDK</div><div>NKDDQITLDE</div><div>FKEAAKSDPS</div><div>IVLLLQCDIQ</div><div>K</div></div>


技术指标	
Our Abpromise guarantee covers the use of ab63276 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	Western blot
形式	Lyophilized

制备和贮存	
稳定性和存储	Shipped at 4°C. Upon delivery aliquot. Store at -80°C. Avoid freeze / thaw cycle. Constituents: 0.242% Tris, 0.116% Sodium chloride
复溶	Add deionized water to prepare a working stock solution of approximately 0.5 mg/mL and let the lyophilized pellet dissolve completely. Product is not sterile! Please filter the product by an appropriate sterile filter before using it in the cell culture.

常规信息

功能	Regulates (in vitro) the inhibition of rhodopsin phosphorylation in a calcium-dependent manner.
组织特异性	Brain and retina. Neuron-specific in the central and peripheral nervous system. Increased in the cerebrospinal fluid of Alzheimer disease patients (at protein level).
序列相似性	Belongs to the recoverin family. Contains 4 EF-hand domains.

图片



LANES:

1. MW marker – 14, 21, 31, 45, 66, 97 kDa
2. reduced and heated sample, 5µg/lane
3. non-reduced and non-heated sample, 5µg/lane

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors