

Recombinant Human VSTM2L protein ab116198

1 图像

描述	
产品名称	重组人VSTM2L蛋白
纯度	> 90 % SDS-PAGE. ab116198 was purified using conventional chromatography.
表达系统	Escherichia coli
Accession	<u>Q96N03</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MGSSHHHHHHSSGLVPRGSHMGSHMTRPAGHAPWDNHVSGHALFTETPHD MTARTGEDVEMACSFRRGSGSPSYSLEIQWYVRSRDWTDKQAWASNQLK ASQQEDAGKEATKISVVKVVGSNISHKLRLSRVKPTDEGSYECRVIDFSD GKARHHKVKAYLRVQPGENSVLHLPEAPPAAPAPPPPKPGKELRKRSVDQ EACSL
预测分子量	23 kDa including tags
氨基酸	25 to 204
标签	His tag N-Terminus

技术指标	
Our Abpromise guarantee covers the use of ab116198 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	Mass Spectrometry SDS-PAGE
质谱法	MALDI-TOF
形式	Liquid

制备和贮存

稳定性和存储

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00


Constituents: 0.002% PMSF, 0.02% DTT, 0.32% Tris HCl, 20% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

常规信息

序列相似性

Contains 1 Ig-like (immunoglobulin-like) domain.

图片



15% SDS-PAGE analysis of 3 µg ab116198.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors