

Recombinant Mouse CXCL13 protein ab50165

描述	
产品名称	重组小鼠CXCL13蛋白
纯度	> 95 % SDS-PAGE. ab50165 purity was assessed also by HPLC. Endotoxin level is less than 0.1 ng per µg (1EU/µg).
内毒素水平	< 0.100 Eu/µg
表达系统	Escherichia coli
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Mouse
序列	<div> <div> ILEAHYTNLK CRCSGVISTV VGLNIIDRIQ VTPPGNGCPK TEVVIWTKMK KVICVNPRAK WLQRLLRHVQ SKSLSSSTPQA PVSKRRAA </div> </div>
氨基酸	22 to 109
描述	重组小鼠BCA1蛋白

技术指标	
Our Abpromise guarantee covers the use of ab50165 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	<div> <div>SDS-PAGE</div> <div>Functional Studies</div> </div>
形式	Lyophilized

制备和贮存	
稳定性和存储	Shipped at 4°C. The lyophilized protein is stable for a few weeks at room temperature. Store at -20°C long term.
复溶	We recommend a quick spin followed by reconstitution in water to a concentration of 0.5 -1 mg/ml. This solution can then be diluted into other aqueous buffers and stored at 4°C for 1 week or -20°C for future use.

常规信息

功能	Chemotactic for B-lymphocytes but not for T-lymphocytes, monocytes and neutrophils. Does not induce calcium release in B-lymphocytes. Binds to BLR1/CXCR5.
组织特异性	Highest levels in liver, followed by spleen, lymph node, appendix and stomach. Low levels in salivary gland, mammary gland and fetal spleen.
序列相似性	Belongs to the intercrine alpha (chemokine CxC) family.
细胞定位	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors