

Recombinant Mouse Eotaxin 2 protein ab50164

描述

产品名称	重组小鼠Eotaxin 2蛋白
纯度	> 95 % SDS-PAGE. ab50164 purity was assessed also by HPLC. Endotoxin level is less than 0.1 ng per µg(1EU/µg).
内毒素水平	< 0.100 Eu/µg
表达系统	Escherichia coli
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Mouse
序列	VTIPSSCCTS FISKKIPENR VVSYQLANGS ICPKAGVIFI TKKGHKICTD PKLLWVQRHI QKLDAKKNQP SKGAKAVRTK FAVQRRRGNS TEV
氨基酸	27 to 119

技术指标

Our **Abpromise guarantee** covers the use of **ab50164** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE Functional Studies
形式	Lyophilized

制备和贮存

稳定性和存储	Shipped at 4°C. The lyophilized protein is stable for a few weeks at room temperature. Store at -20°C long term.
复溶	Reconstitute to 1mg/ml in distilled water.

常规信息

功能	Chemotactic for resting T-lymphocytes, and eosinophils. Has lower chemotactic activity for
----	--

neutrophils but none for monocytes and activated lymphocytes. Is a strong suppressor of colony formation by a multipotential hematopoietic progenitor cell line. Binds to CCR3.

组织特异性	Activated monocytes and activated T lymphocytes.
序列相似性	Belongs to the intercrine beta (chemokine CC) family.
翻译后修饰	N-glycosylated.
细胞定位	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors