

Recombinant mouse MCP3 protein ab9912

描述

产品名称	重组小鼠MCP3蛋白
生物活性	Determined by its ability to chemoattract BaF3 cells expressing human CCR2 receptor using a concentration range of 300-500 ng/mL.
纯度	>= 98 % SDS-PAGE. = 98% by SDS-PAGE gel and HPLC analyses.s.
内毒素水平	< 1.000 Eu/μg
表达系统	Escherichia coli
Accession	<u>Q03366</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Mouse
序列	QPDGPNASTCCYVKKQKIPKRNLKSYRRITSSRCPWEAVIFK TKKGMEVC AEAHQKWVEEAIAYLDMKTPTPKP
预测分子量	11 kDa
氨基酸	24 to 97

技术指标

Our **Abpromise guarantee** covers the use of **ab9912** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Functional Studies
	SDS-PAGE
形式	Lyophilized

制备和贮存

稳定性和存储	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. Sterile filtered through a 0.2 micron filter. Lyophilized with no additives. This product is an active protein and may elicit a biological response in vivo, handle with caution.
--------	--

复溶

For lot specific reconstitution information please contact our Scientific Support Team.

常规信息

功能	Chemotactic factor that attracts monocytes and eosinophils, but not neutrophils. Augments monocyte anti-tumor activity. Also induces the release of gelatinase B. This protein can bind heparin. Binds to CCR1, CCR2 and CCR3.
序列相似性	Belongs to the intercrine beta (chemokine CC) family.
翻译后修饰	O-glycosylated.
细胞定位	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors