

# Recombinant Mouse Pancreatic Lipase/PTL protein (His tag) ab226285

### 1 图像

#### 描述

<b>产品名称</b>	重组小鼠Pancreatic Lipase/PTL蛋白(His tag)
<b>纯度</b>	> 90 % SDS-PAGE.
<b>表达系统</b>	Escherichia coli
<b>Accession</b>	<u><a href="#">Q6P8U6</a></u>
<b>蛋白长度</b>	Full length protein
<b>无动物成分</b>	No
<b>性质</b>	Recombinant
<b>种属</b>	Mouse
<b>序列</b>	<pre> REVCFDKLGCFSDDAPWSGTLDRPLKALPWSPAQINTRFLLY TNENPDNY QLITSDASNIRNSNFRITNRKTRIIIHGFIDKGEENWLSDMCK NMF RVESV NCICVDWKGGSRTTYTQATQNV RVVGAEVALLVNVLQSDLGY SLNNVHLI GHS LGSHIAGEAGKRTFGAIGRITGLDPAEPYFQGTPEEVRL DPTDAQFV DAIHTDAGPIIPNLGFGMSQTVGHLDFFPNGGIEMPGCQKNI LSQIVDID GIWEGTRNFAACNHLRSYKFYTDSIVNPTGFAGFSCSSYSLF TANKCFPC GSGGCPQMGHYADRYPGKTSRLYQTFYLN TGDKSNFARWRYQ VTVTLSGQ KVTGHILVSLFGNGGNSKQYEVFKGSLQPGTSHVNEFDSVDV VGDLQKVK FIWYNNVINPTLPKVGASRITVERNDGRVFNFCSETVREDV LLTLSPC </pre>
<b>预测分子量</b>	66 kDa including tags
<b>氨基酸</b>	17 to 465
<b>标签</b>	His tag N-Terminus
<b>额外的序列信息</b>	This product is the mature full length protein from aa 17 to 465 with 6xHis-SUMO tag at the N-

terminus. The signal peptide is not included.

## 技术指标

---

Our **Abpromise guarantee** covers the use of **ab226285** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE
形式	Liquid
补充说明	Previously labelled as Pancreatic Lipase.

## 制备和贮存

---

稳定性和存储	Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. pH: 7.2 Constituents: 50% Glycerol (glycerin, glycerine), Tris buffer
--------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------

## 常规信息


---

序列相似性	Belongs to the AB hydrolase superfamily. Lipase family. Contains 1 PLAT domain.
细胞定位	Secreted.

---

## 图片

---



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) analysis of ab226285 with 5% enrichment gel and 15% separation gel.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- 
- Replacement or refund for products not performing as stated on the datasheet
  - Valid for 12 months from date of delivery
  - Response to your inquiry within 24 hours
  
  - We provide support in Chinese, English, French, German, Japanese and Spanish
  - Extensive multi-media technical resources to help you
  - We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors