

Recombinant mouse PDGF B protein ab78593

2 图像

描述

产品名称	重组小鼠PDGF B蛋白
纯度	= 11 % SDS-PAGE. Purity Greater than 98% as determined by: Reducing and non-reducing SDS-PAGE. Analytical HPLC. Endotoxin level as measured by LAL is <0.01ng/ug or <0.1EU/ug. This product was sterile filtered.
表达系统	Escherichia coli
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Mouse
序列	MSLGSLAAAE PAVIAECKTR TEVFQISRNL IDRTNANFLV WPPCVEVQRC SGCCNNRNVQ CRASQVQMRP VQVRKIEIVR KKPIFKKATV TLEDHLACKC ETIVTPRPVT

技术指标

Our **Abpromise guarantee** covers the use of **ab78593** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	SDS-PAGE
形式	Lyophilized
补充说明	It is recommended to add a carrier protein (0.1% HSA or BSA) for long term storage.

制备和贮存

稳定性和存储	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. Constituent: 0.294% Sodium citrate This product is an active protein and may elicit a biological response in vivo, handle with caution.
复溶	Reconstitute the lyophilized product with sterile water at a concentration of 0.1 - 0.5 mg's/ml, which can be further diluted into other aqueous solutions.


常规信息

图片



SDS-PAGE - Recombinant mouse PDGF BB protein (ab78593)

SDS PAGE analysis of ab78593 under non-reducing (-) and reducing (+) conditions. Stained with Coomassie Blue.



Functional Studies - Recombinant mouse PDGF BB protein (ab78593)

ab78593 used in Functional Studies.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors