

# Recombinant Mouse SF20 protein ab54228

### 描述

---

产品名称	重组小鼠SF20蛋白
纯度	> 95 % SDS-PAGE. Purity: > 98% by SDS-PAGE
表达系统	Escherichia coli
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Mouse

### 技术指标

---

Our **Abpromise guarantee** covers the use of **ab54228** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE
形式	Lyophilized

### 制备和贮存

---

稳定性和存储	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
复溶	Reconstitute in water to a concentration of 0.1-1.0 mg/ml. This solution can be diluted in water or other buffer solutions.

### 常规信息

---

组织特异性	Expressed in synovial tissue. Found in synovial fluid of patients with arthropaties.
序列相似性	Belongs to the UPF0556 family.
细胞定位	Secreted. Endoplasmic reticulum-Golgi intermediate compartment.

---

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors