

Recombinant Mouse Vaspin protein ab109036

描述

产品名称	重组小鼠Vaspin蛋白
纯度	> 90 % SDS-PAGE. ab109036 is 0.2µm filtered.
内毒素水平	< 0.100 Eu/µg
表达系统	HEK 293 cells
Accession	<u>Q7TMF5</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Mouse
预测分子量	50 kDa including tags
氨基酸	21 to 413
标签	DDDDK tag N-Terminus

技术指标

Our **Abpromise guarantee** covers the use of **ab109036** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE
形式	Liquid
补充说明	This product was previously labelled as Vaspin, SERPINA12

制备和贮存

稳定性和存储	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. pH: 7.20 Constituent: PBS
--------	---

常规信息

功能	May modulates insulin action conceivably only in the presence of its yet undefined target proteases in white adipose tissues.
组织特异性	Expressed in visceral adipose tissues.
序列相似性	Belongs to the serpin family.
细胞定位	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors