

# Recombinant mouse Visfatin protein ab62617

### 描述

---

产品名称	重组小鼠Visfatin蛋白
生物活性	Biological activity: Mimics insulin effects in differentiating 3T3-L1 cells.
纯度	> 95 % SDS-PAGE. Purity >95% as determined by SDS-PAGE and HPLC ENDOTOXIN: <0.1 EU/μg protein (LAL test)
表达系统	Escherichia coli
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Mouse

### 技术指标

---

Our **Abpromise guarantee** covers the use of **ab62617** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE Functional Studies
形式	Lyophilized
补充说明	This product is manufactured by BioVision, an Abcam company and was previously called 4908 Visfatin, mouse recombinant. 4908-10 is the same size as the 10 μg size of ab62617.

### 制备和贮存

---

稳定性和存储	Shipped at 4°C. Upon delivery aliquot. Store at -80°C. Avoid freeze / thaw cycle. Constituent: PBS This product is an active protein and may elicit a biological response in vivo, handle with caution.
复溶	Reconstitute in PBS to 1mg/ml

### 常规信息

---

<b>功能</b>	Catalyzes the condensation of nicotinamide with 5-phosphoribosyl-1-pyrophosphate to yield nicotinamide mononucleotide, an intermediate in the biosynthesis of NAD. It is the rate limiting component in the mammalian NAD biosynthesis pathway.
<b>组织特异性</b>	Expressed in large amounts in bone marrow, liver tissue, and muscle. Also present in heart, placenta, lung, and kidney tissues.
<b>通路</b>	Cofactor biosynthesis; NAD(+) biosynthesis; nicotinamide D-ribonucleotide from 5-phospho-alpha-D-ribose 1-diphosphate and nicotinamide: step 1/1.
<b>序列相似性</b>	Belongs to the NAPRTase family.
<b>细胞定位</b>	Cytoplasm.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors