


Recombinant Polyomavirus strain RA Major Capsid VP1 protein ab74571

1 图像

描述	
产品名称	重组Polyomavirus strain RA Major Capsid VP1蛋白
纯度	> 75 % SDS-PAGE. Purified by ultracentrifugation
表达系统	Saccharomyces cerevisiae
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
技术指标	
Our Abpromise guarantee covers the use of ab74571 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	ELISA Western blot SDS-PAGE
形式	Lyophilized
制备和贮存	
稳定性和存储	Shipped at 4°C. Store at +4°C. Constituent: PBS
复溶	Reconstitute with deionized H ₂ O. After reconstitution store at 4°C.
常规信息	
相关性	The polyomavirus capsid is built of 72 capsomeres that are all pentamers of the major capsid protein VP. The 12 peavalent and 60 hexavalent capsomeres are arranged in the positions of the 5- and 6-coordinated vertices of a T = 7d icosahedral surface lattice.
细胞定位	Virion. Nucleus

图片



SDS-PAGE - Recombinant Polyomavirus strain RA
Major Capsid VP1 protein (ab74571)

SDS-PAGE showing ab74571 at approximately 43kDa
(4.4µg/lane).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors