

Recombinant Rat CINC1 protein ab54466

描述

产品名称	重组大鼠CINC1蛋白
纯度	> 95 % SDS-PAGE. Purity: > 98% (SDS-PAGE; HPLC).
表达系统	Escherichia coli
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Rat

技术指标

Our **Abpromise guarantee** covers the use of **ab54466** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE Functional Studies
形式	Lyophilized

制备和贮存

稳定性和存储	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
复溶	The lyophilized rat CINC1 is soluble in water and most aqueous buffers. The lyophilized powder should be reconstituted in water to a concentration of 100µg/ml. This solution can be diluted into water or other buffered solutions.

常规信息

相关性	Rat CINC1 (Cytokine Induced Neutrophil Chemoattractant 1) is a member of the CXC subfamily of chemokines which also includes three other rat CXC chemokines (CINC2 alpha, CINC2 beta, CINC3/MIP2). CINC1 shares 63 % to 67 % amino acid sequence identity with the other rat CINC1s. Based on their sequence homology, it has been suggested that CINC1s are the rat counterpart of the human GROs. Rat CINC1s, potent neutrophil attractants and activators, have an important role in the infiltration of neutrophils into inflammatory sites in rats.
-----	--

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors