

Recombinant Rat ErbB2 / HER2 protein ab191967

1 References

描述	
产品名称	重组大鼠ErbB2 / HER2蛋白
纯度	> 95 % SDS-PAGE. Purity greater than 95% as determined by SEC-HPLC and reducing SDS-PAGE.
内毒素水平	< 1.000 Eu/μg
表达系统	Escherichia coli
Accession	<b><u>P06494</u></b>
蛋白长度	Protein fragment
无动物成分	No
性质	Recombinant
种属	Rat
序列	MANASLSFLQDIQEVQGYMLIAHNQVKRVPLQRLRIVRGTQL FEDKYALA VLDNRDPQDNVAASTPGRTPEGLRELQLRSLTEILKGGVLIR GNPQLCYQ DMVLWKDVFRKNNQLAPVDIDTNRSRACPPCAPACKDNHCWG ESPEDCQI LTGTICTSGCARCKGRLPTDCCHEQCAAGCTGPKHSDCLACL HFNHSGIC ELHCPALVTYNTDTFESMHNPEGRYTFGASCVTTCPYNYLST EVGSCTLV CPPNNQEVHHHHHH
预测分子量	29 kDa including tags
氨基酸	67 to 323
标签	His tag C-Terminus

技术指标	
Our <b>Abpromise guarantee</b> covers the use of <b>ab191967</b> in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	SDS-PAGE HPLC
形式	Lyophilized

## 制备和贮存

---

### 稳定性和存储

Shipped at 4°C. The lyophilized protein is stable for a few weeks at room temperature. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 7.4

Constituents: 0.87% Sodium chloride, 24% Urea, 0.24% Tris

Lyophilized from a 0.2 µM filtered solution.

## 常规信息

---

### 功能

Protein tyrosine kinase that is part of several cell surface receptor complexes, but that apparently needs a coreceptor for ligand binding. Essential component of a neuregulin-receptor complex, although neuregulins do not interact with it alone. GP30 is a potential ligand for this receptor. Regulates outgrowth and stabilization of peripheral microtubules (MTs). Upon ERBB2 activation, the MEMO1-RHOA-DIAPH1 signaling pathway elicits the phosphorylation and thus the inhibition of GSK3B at cell membrane. This prevents the phosphorylation of APC and CLASP2, allowing its association with the cell membrane. In turn, membrane-bound APC allows the localization of MACF1 to the cell membrane, which is required for microtubule capture and stabilization. In the nucleus is involved in transcriptional regulation. Associates with the 5'-TCAAATTC-3' sequence in the PTGS2/COX-2 promoter and activates its transcription. Implicated in transcriptional activation of CDKN1A; the function involves STAT3 and SRC. Involved in the transcription of rRNA genes by RNA Pol I and enhances protein synthesis and cell growth.

### 组织特异性

Expressed in a variety of tumor tissues including primary breast tumors and tumors from small bowel, esophagus, kidney and mouth.

### 疾病相关

Hereditary diffuse gastric cancer

Glioma

Ovarian cancer

Lung cancer

Gastric cancer

Chromosomal aberrations involving ERBB2 may be a cause gastric cancer. Deletions within 17q12 region producing fusion transcripts with CDK12, leading to CDK12-ERBB2 fusion leading to truncated CDK12 protein not in-frame with ERBB2.

### 序列相似性

Belongs to the protein kinase superfamily. Tyr protein kinase family. EGF receptor subfamily.

Contains 1 protein kinase domain.

### 翻译后修饰

Autophosphorylated. Autophosphorylation occurs in trans, i.e. one subunit of the dimeric receptor phosphorylates tyrosine residues on the other subunit (Probable). Ligand-binding increases phosphorylation on tyrosine residues (PubMed:27134172). Signaling via SEMA4C promotes phosphorylation at Tyr-1248 (PubMed:17554007). Dephosphorylated by PTPN12 (PubMed:27134172).

### 细胞定位

Cytoplasm. Nucleus and Cell membrane. Cytoplasm, perinuclear region. Nucleus. Translocation to the nucleus requires endocytosis, probably endosomal sorting and is mediated by importin beta-1/KPNB1.

---

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors