

Recombinant Rat Prolactin/PRL protein ab54140

描述

| | |
|-------|---|
| 产品名称 | 重组大鼠Prolactin/PRL蛋白 |
| 纯度 | > 95 % SDS-PAGE. Purity : Greater than 98% by SDS-PAGE gel and HPLC analyses. |
| 内毒素水平 | < 0.100 Eu/μg |
| 表达系统 | Escherichia coli |
| 蛋白长度 | Full length protein |
| 无动物成分 | No |
| 性质 | Recombinant |
| 种属 | Rat |
| 序列 | MLPVCSGGDC QTPLPELFDR VVMLSHYIHT LYTDMFIEFD KQYVQDREFI AKAINDCPTS SLATPEDKEQ AQKVPPEVLL NLILSLVHSW NDPLFQLITG LGGIHEAPDA IISRAKEIEE QNKRLLEGIE KIISQAYPEA KGNEIYLVWS QLPSLQGVDE ESKDLAFYNN IRCLRRDSHK VDNYLKFLRC QIVHKNNC |

技术指标

Our **Abpromise guarantee** covers the use of **ab54140** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| | |
|------|--|
| 应用 | SDS-PAGE Functional Studies |
| 形式 | Lyophilized |
| 补充说明 | Reconstituted Prolactin is stable for at least 3 months. |

制备和贮存

| | |
|--------|--|
| 稳定性和存储 | Shipped at 4°C. Upon delivery aliquot. Upon reconstitution add a carrier protein (0.1% BSA). Store at -20°C. Avoid freeze / thaw cycle. For long term storage it is recommended to add a carrier protein on reconstitution (0.1% HSA or BSA). Please see notes section. |
| 复溶 | For lot specific reconstitution information please contact our Scientific Support Team. |

常规信息

| | |
|------|--|
| 相关性 | Various hormones are secreted from the anterior pituitary during development and growth. Prolactin is a growth factor secreted by the anterior pituitary that is necessary for the proliferation and differentiation of the mammary glands. |
| 细胞定位 | Secreted |

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors