

450 nm Stop Solution for TMB Substrate ab171529

17 References

概述

产品名称	450 nm Stop Solution for TMB底物
经测试应用	适用于: ELISA
常规说明	<p>Abcam's 450 nm Stop Solution for TMB Substrate is a ready-to-use solution for colorimetric HRP-based ELISA detection. The acidic formulation is a proprietary solution used to stop TMB ELISA Substrate reactions (ab171522, ab171523, ab171524, ab171525, ab171526 and ab171527) and changes the chromogen from blue (620 nm) to yellow (450 nm) in color.</p> <p>Product Use:</p> <p>Equal volumes of TMB ELISA Substrate and 450 nm Stop Solution for TMB Substrate should be used.</p>

性能

形式	Liquid
存放说明	Shipped at Room Temperature. Store at Room Temperature.
存储溶液	<p>pH: 1</p> <p>Constituent: 5% Methanesulphonic acid</p>

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab171529于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		Use at an assay dependent concentration. Equal volumes of TMB ELISA Substrate and 450 nm Stop Solution for TMB Substrate should be used.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors