

10kD Spin Column ab93349

★★★★★ [1 Abreviews](#) [27 References](#) [1 图像](#)

概述

产品名称	10kD Spin Column
产品概述	<p>The 10kD Spin Columns (ab93349) are disposable ultrafiltration devices for the concentration or deproteinization of biological samples.</p> <p>They can be used in either a swing bucket or fixed angle rotors which accept 2.0 mL centrifuge tube at maximum speed 10000 x <i>g</i> (5-10 mins). The vertical membrane design and thin channel filtration chamber minimizes membrane fouling and provides high speed concentration, even laden solutions. It is specifically designed with low internal surface and membrane area in order to achieve superior recovery from dilute solution.</p> <p>Maximum initial sample volume is 500µL. The dead stop volume is approximately 20-25µL.</p>
说明	<p>This product is manufactured by BioVision, an Abcam company and was previously called 1997 10K Spin Column. 1997-25 is the same size as the 25 unit size of ab93349.</p> <p>If working with cell or tissue lysates, we generally recommend the <a href="#">TCA Sample Deproteinization Kit</a> or using a <a href="#">PCA/KOH Deproteinization protocol</a>. If it is critical to use a spin column with a cell or tissue lysate then, to avoid blocking the filter, first centrifuge the lysate samples at ~14,000 rpm for 10-20 mins to pellet debris, and then use the resulting supernatant with the spin column.</p>
经测试应用	适用于: Purification

性能

存放说明 Store at room temperature. Please refer to protocols.

组件	25 units	25 units
Concentrator cup: Polycarbonate	25 units	25 units
Filtrate vessel: Polypropylene	25 units	25 units
Membrane: 10,000 MWCO Polyethersulfone	25 units	25 units

## 应用

### The Abpromise guarantee

**Abpromise™**承诺保证使用ab93349于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Purification		Use at an assay dependent concentration.

## 图片



10kD Spin Column (ab93349)

Representative 10kD Spin Column (ab93349).

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors