

Donkey Anti-Chicken IgY H&L (FITC) ab63507

★★★★★ [2 Abreviews](#) [15 References](#) [2 图像](#)

概述

产品名称	驴抗鸡IgY H&L (FITC)
宿主	Donkey
靶标种属	Chicken
经测试应用	适用于: IHC-Fr
免疫原	Purified chicken IgY isolated from normal chicken egg yolks.
偶联物	FITC. Ex: 493nm, Em: 528nm

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.20 Preservative: 0.01% Thimerosal (merthiolate) Constituents: PBS, 50% Glycerol (glycerin, glycerine)
纯度	Protein L purified
克隆	多克隆
同种型	IgG


应用

The Abpromise guarantee **Abpromise™** 承诺保证使用ab63507于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-Fr	★★★★★ (1)	Use at an assay dependent dilution.

图片




Immunohistochemistry (Frozen sections) - Donkey Anti-Chicken IgY H&L (FITC) (ab63507)

This image is courtesy of an anonymous Abreview.

ab13970 staining GFP in E10.5 mouse embryo sections by Immunohistochemistry (IHC-Fr - frozen sections).

Tissue was fixed with paraformaldehyde, permeabilized with PBS/0.5% Triton-X100 and blocked with protein block for 30 minutes at 25°C. Samples were incubated with primary antibody (1/500) for 3 hours at 25°C. FITC-conjugated donkey anti-chicken IgY H&L (ab63507) (1/500) was used as the secondary antibody.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Donkey Anti-Chicken IgY H&L (FITC) (ab63507)

ab63507 was used at dilution 1/100 with the primary antibody **ab13970** in IHC - Free Floating. See the review on **ab13970**.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors