

Donkey Anti-Goat IgG H&L (HRP polymer) ab214881

★★★★★ [1 Abreviews](#) [1 图像](#)

概述

产品名称	驴抗山羊IgG H&L (HRP polymer)
宿主	Donkey
靶标种属	Goat
特异性	Based on IEP, ab214881 reacts with heavy (gamma) chains on goat IgG, light chains on all goat immunoglobulins. No cross-reactivity is observed to non-immunoglobulin goat serum proteins, IgG from human, mouse, rabbit or rat.
经测试应用	适用于: Immunoelectrophoresis, IHC-P
免疫原	Full length protein corresponding to Goat IgG. (Purified Goat IgG, whole molecule).
偶联物	HRP polymer

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.1% Proclin 300 Constituents: 1% BSA, 98% PBS
纯度	Affinity purified
纯化说明	ab214881 was affinity purified using solid phase Goat IgG. Purity >95% by SDS-PAGE.
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee


Abpromise™ 承诺保证使用ab214881于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明

应用	Ab评论	说明
Immunoelectrophoresis		Use at an assay dependent concentration. This product is provided in a pre-diluted, ready-to-use format. Incubate 150 µl per slide for 1 hour.
IHC-P	★★★★★ (1)	Use at an assay dependent concentration. This product is provided in a pre-diluted, ready-to-use format. Incubate 150 µl per slide for 1 hour.

图片



Human tonsil tissue stained with 1 µg/ml goat anti-human IgG (H&L) detected with ab214881 Donkey Anti-Goat H&L (HRP polymer).

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Donkey Anti-Goat IgG H&L (HRP polymer) (ab214881)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors