

# Donkey Anti-Mouse IgG H&L (Alkaline Phosphatase) ab97027

## 5 References

### 概述

|       |                                                                                                                                                                                                                     |
|-------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 产品名称  | 驴抗小鼠IgG H&L (Alkaline Phosphatase)                                                                                                                                                                                  |
| 宿主    | Donkey                                                                                                                                                                                                              |
| 靶标种属  | Mouse                                                                                                                                                                                                               |
| 特异性   | By immunoelectrophoresis and ELISA this antibody reacts specifically with mouse IgG and with light chains common to other mouse immunoglobulins. No antibody was detected against nonimmunoglobulin serum proteins. |
| 经测试应用 | 适用于: ICC, IHC-P, ELISA, WB                                                                                                                                                                                          |
| 偶联物   | Alkaline Phosphatase                                                                                                                                                                                                |

### 性能

|      |                                                                                                                                                               |
|------|---------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 形式   | Liquid                                                                                                                                                        |
| 存放说明 | Shipped at 4°C. Store at +4°C.                                                                                                                                |
| 存储溶液 | pH: 6.8<br>Preservative: 0.09% Sodium azide<br>Constituents: 0.2% BSA, 0.58% Sodium chloride, 0.0095% Magnesium chloride, 0.00136% Zinc chloride, 1.19% HEPES |
| 纯度   | Immunogen affinity purified                                                                                                                                   |
| 纯化说明 | This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to Alkaline Phosphatase (AP).                     |
| 共轭说明 | Molar enzyme: antibody protein ratio is 1:1                                                                                                                   |
| 克隆   | 多克隆                                                                                                                                                           |
| 同种型  | IgG                                                                                                                                                           |

### 应用

**The Abpromise guarantee**      **Abpromise™** 承诺保证使用ab97027于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

| 应用    | Ab评论 | 说明                                   |
|-------|------|--------------------------------------|
| ICC   |      | Use at an assay dependent dilution.  |
| IHC-P |      | 1/100 - 1/1000.                      |
| ELISA |      | 1/5000 - 1/50000. (Primary)          |
| WB    |      | 1/5000 - 1/50000. (Chemiluminescent) |

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors