

# Donkey Anti-Mouse IgG H&L (DyLight® 488) ab96875

[18 References](#)   [2 图像](#)

### 概述

产品名称	驴抗小鼠IgG H&L (DyLight® 488)
宿主	Donkey
靶标种属	Mouse
特异性	By immunoelectrophoresis and ELISA this antibody reacts specifically with Mouse IgG and with light chains common to other Mouse immunoglobulins. No antibody was detected against non immunoglobulin serum proteins.
经测试应用	<b>适用于:</b> ICC/IF, IHC-P, Flow Cyt
免疫原	Other Immunogen Type corresponding to Mouse IgG.
偶联物	DyLight® 488. Ex: 493nm, Em: 518nm

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: 0.2% BSA, PBS
纯度	Immunogen affinity purified
纯化说明	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to DyLight® 488.
克隆	多克隆
同种型	IgG


### 应用

**The Abpromise guarantee**      **Abpromise™** 承诺保证使用ab96875于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		1/50 - 1/500.
IHC-P		1/50 - 1/500.
Flow Cyt		1/50 - 1/200.

## 图片




Immunocytochemistry/ Immunofluorescence -  
Donkey Anti-Mouse IgG H&L (DyLight® 488)  
(ab96875)

This image is courtesy of an Abreview submitted by Charles Pallangyo.

**ab7817** staining alpha smooth muscle actin in mouse primary colon myofibroblasts by ICC/IF

(Immunocytochemistry/immunofluorescence). Cells were fixed with acetone and blocked with 5% BSA for 30 minutes at 25°C. Samples were incubated with primary antibody (1/100 in PBS + 5% BSA) for 2 hour at 25°C. A DyLight® 488-conjugated donkey anti-mouse IgG H&L (ab96875) (1/1000) was used as the secondary antibody. Vimentin appears in red.



Donkey Anti-Mouse IgG H&L (DyLight® 488)  
(ab96875)

Emission spectra of DyLight® fluorochromes available in our catalog.

Line colors represent the approximate visible colors of the wavelength maxima.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors