

# Donkey Anti-Mouse IgG H&L (HRP) preadsorbed ab98799

## 8 References

### 概述

|        |  |                              |
|--------|--|------------------------------|
| 产品名称   | 驴抗小鼠IgG H&L (HRP)预吸附二抗                               |                              |
| 宿主     | Donkey   |                              |
| 靶标种属   | Mouse  |                              |
| 经测试应用  | 适用于: IHC-P, ELISA, WB, ICC                           |                              |
| 最小交叉反应 | Chicken, Cow, Goat, Horse, Human, Rabbit, Rat, Sheep | <a href="#">more details</a> |
| 偶联物    | HRP  |                              |

### 性能

|      |   |
|------|---|
| 形式   | Liquid  |
| 存放说明 | Shipped at 4°C. Store at +4°C.  |
| 存储溶液 | pH: 6.8<br>Constituents: 0.2% BSA, 0.05% CMIT/MIT based preservative  |
| 纯度   | Immunogen affinity purified   |
| 纯化说明 | This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to Horse Radish Peroxidase. |
| 克隆   | 多克隆   |
| 同种型  | IgG   |

### 应用

**The Abpromise guarantee**      **Abpromise™** 承诺保证使用ab98799于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

| 应用    | Ab评论 | 说明                                    |
|-------|------|---------------------------------------|
| IHC-P |      | 1/200 - 1/5000.                       |
| ELISA |      | 1/10000 - 1/100000. Primary dilution. |

| 应用  | Ab评论 | 说明   |
|-----|------|--|
| WB  |      | 1/5000 - 1/50000. Colorimetric: Use 1/5,000-1/30,000 dilution; Chemiluminescent: Use 1/10,000-1/50,000 dilution. |
| ICC |      | Use at an assay dependent dilution.  |

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors