

Donkey Anti-Mouse IgG H&L (PE) preadsorbed ab7003

9 References 1 图像

概述	
产品名称	驴抗小鼠IgG H&L (PE)预吸附二抗
宿主	Donkey
靶标种属	Mouse
经测试应用	适用于: IHC-P, IHC-Fr, Immunomicroscopy, Flow Cyt, ICC/IF
最小交叉反应	Chicken, Cow, Goat, Guinea pig, Hamster, Horse, Human, Rabbit, Rat, Sheep <a href="#">more details</a>
免疫原	Mouse IgG whole molecule
偶联物	PE. Ex: 488nm, Em: 575nm
性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 6.50 Preservative: 0.01% Sodium azide Constituents: 0.42% Tripotassium orthophosphate, 0.87% Sodium chloride, 1% BSA
纯度	Immunogen affinity purified
纯化说明	Prepared from monospecific antiserum by immunoaffinity chromatography using Mouse IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities, pepsin digestion and chromatographic separation
共轭说明	R-Phycoerythrin (R-PE) Molecular Weight 240,000 daltons) Absorption Wavelength: 488 nm Emission Wavelength: 575 nm OD <sub>566nm</sub> /OD <sub>280nm</sub> : 2.8
克隆	多克隆
同种型	IgG
常规说明	Can be used in immunomicroscopy and flow cytometry or FACS analysis as well as other antibody based fluorescent assays requiring extremely low background levels, absence of F(c) mediated binding, lot-to-lot consistency, high titer and specificity.

## 应用


### The Abpromise guarantee

**Abpromise™**承诺保证使用ab7003于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P		Use at an assay dependent dilution.
IHC-Fr		Use at an assay dependent dilution.
Immunomicroscopy		Use at an assay dependent dilution.
Flow Cyt		1/100 - 1/250.
ICC/IF		1/100 - 1/250.

## 图片



ab7003 was used with the primary antibody **ab33315** in Flow Cyt.

See Abreview on **ab33315**.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors