

Donkey Anti-Rabbit IgG H&L (12nm Gold) preadsorbed  
ab105295

1 References

概述	
产品名称	驴抗兔IgG H&L (12nm Gold)预吸附二抗
宿主	Donkey
靶标种属	Rabbit
特异性	Recognizes rabbit IgG, H&L. X-Adsorbed and exhibits minimal crossreactivity to bovine, chicken, goat, guinea pig, Syrian hamster, equine, human, mouse, rat and sheep serum proteins.
经测试应用	适用于: Electron Microscopy, ELISA, WB
最小交叉反应	Chicken, Cow, Goat, Guinea pig, Horse, Human, Mouse, Rat, Sheep, Syrian hamster
<a href="#">more details</a>	
偶联物	Gold 12nm

性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 8.50 Preservative: 0.05% Sodium azide Constituents: 0.164% Sodium phosphate, 0.381% Sodium borate, 1.5% BSA, 0.87% Sodium chloride
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用	
The Abpromise guarantee	<a href="#">Abpromise™</a> 承诺保证使用ab105295于以下的经测试应用
“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。	

应用	Ab评论	说明
Electron Microscopy		1/20 - 1/40.
ELISA		Use at an assay dependent dilution.
WB		Use at an assay dependent dilution.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors