

Donkey Anti-Rabbit IgG H&L (Texas Red ®) preadsorbed ab7081

5 References

概述

| | |
|--------|---|
| 产品名称 | 驴抗兔IgG H&L (Texas Red ®)预吸附二抗 |
| 宿主 | Donkey |
| 靶标种属 | Rabbit |
| 经测试应用 | 适用于: ICC/IF, Immunomicroscopy, Flow Cyt, IHC-P, IHC-Fr, ELISA |
| 最小交叉反应 | Chicken, Cow, Goat, Guinea pig, Hamster, Horse, Human, Mouse, Rat, Sheep more details |
| 免疫原 | Full length native Rabbit IgG (purified). |
| 偶联物 | Texas Red ®. Ex: 596nm, Em: 620nm |

性能

| | |
|------|--|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. |
| 存储溶液 | pH: 6.50 Preservative: 0.01% Sodium azide Constituents: 0.42% Tripotassium orthophosphate, 0.87% Sodium chloride, 1% BSA |
| 纯度 | Affinity purified |
| 纯化说明 | Texas Red Sulfonyl Chloride (Molecular Weight 625 daltons) Absorption Wavelength: 596 nm Emission Wavelength: 620 nm Fluorochrome/Protein Ratio: 1.7 moles Texas Red per mole of Goat IgG |
| 共轭说明 | This product was prepared from monospecific antiserum by immunoaffinity chromatography using Rabbit IgG coupled to agarose |
| 克隆 | 多克隆 |

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab7081于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

| 应用 | Ab评论 | 说明 |
|------------------|------|----|
| ICC/IF | | |
| Immunomicroscopy | | |
| Flow Cyt | | |
| IHC-P | | |
| IHC-Fr | | |
| ELISA | | |

应用说明

Flow Cyt: Use at an assay dependant concentration.

IHC-P: Use at an assay dependant concentration.

IHC-Fr: Use at an assay dependant concentration.

IM: Use at an assay dependant concentration.

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors