

Donkey anti-Rat IgG H&L (PE) ab7008

★★★★★ 1 Abreviews 2 图像

概述	
产品名称	驴抗大鼠IgG H&L (PE)
宿主	Donkey
靶标种属	Rat
经测试应用	适用于: ICC/IF, Immunomicroscopy, Flow Cyt, IHC-P, IHC-Fr
最小交叉反应	Chicken, Cow, Goat, Guinea pig, Hamster, Horse, Human, Mouse, Rabbit, Sheep more details
免疫原	Rat IgG whole molecule
偶联物	PE. Ex: 488nm, Em: 575nm

性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 6.50 Preservative: 0.01% Sodium azide Constituents: 0.878% Sodium chloride, 1% BSA, 0.424% Tripotassium orthophosphate
纯度	Immunogen affinity purified
纯化说明	Prepared from monospecific antiserum by immunoaffinity chromatography using rat IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities, pepsin digestion and chromatographic separation
共轭说明	R-Phycoerythrin (R-PE) Molecular Weight 240,000 daltons) Absorption Wavelength: 488 nm Emission Wavelength: 575 nm OD _{566nm} /OD _{280nm} : 2.8
克隆	多克隆
常规说明	Can be used in immunomicroscopy and flow cytometry or FACS analysis as well as other antibody based fluorescent assays requiring extremely low background levels, absence of F(c) mediated binding, lot-to-lot consistency, high titer and specificity.

应用


The Abpromise guarantee

Abpromise™ 承诺保证使用ab7008于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF	★★★★★ (1)	1/100 - 1/250.
Immunomicroscopy		1/50 - 1/200.
Flow Cyt		1/100 - 1/250. The maximum amount of reagent required to stain 1×10^6 cells in flow cytometry is approximately 1.0 ug of antibody conjugate. Lesser amounts of reagent may be sufficient for staining.
IHC-P		1/50 - 1/200.
IHC-Fr		1/50 - 1/200.

图片




Immunocytochemistry/Immunofluorescence analysis of lung tissue using Donkey anti-rabbit IgG H&L (PE) secondary antibody (ab7008).

Immunocytochemistry/ Immunofluorescence -

Donkey anti-Rat IgG H&L (PE) (ab7008)

This image is courtesy of an anonymous Abreview



ab7008 was used with the primary antibody **ab19765** in Flow Cyt.
See Abreview on **ab19765**.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors