

Donkey Anti-Sheep IgG H&L (Alkaline Phosphatase) ab6901

4 References

概述

产品名称	驴抗羊IgG H&L (Alkaline Phosphatase)
宿主	Donkey
靶标种属	Sheep
经测试应用	适用于: ICC/IF, IHC-Fr, IHC-P, WB, ELISA, Immunomicroscopy, Dot blot
免疫原	Full length native Sheep IgG (purified).
偶联物	Alkaline Phosphatase

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.01% Sodium azide Constituents: 0.00136% Zinc chloride, 0.0095% Magnesium chloride, 0.79% Tris HCl, 50% Glycerol (glycerin, glycerine), 0.87% Sodium chloride, 1% BSA
纯度	Affinity purified
纯化说明	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Sheep IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities
共轭说明	Alkaline Phosphatase (Calf Intestine) (Molecular Weight 140,000 daltons)
克隆	多克隆

应用

The Abpromise guarantee Abpromise™承诺保证使用ab6901于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		

应用	Ab评论	说明
IHC-Fr		
IHC-P		
WB		
ELISA		
Immunomicroscopy		
Dot blot		

应用说明

Dot: Use at an assay dependant dilution.

ELISA: 1/3000 - 1/14000.

IHC-P: Use at an assay dependant dilution.

IHC-Fr: Use at an assay dependant dilution.

IM: Use at an assay dependant dilution.

WB: Use at an assay dependant dilution.

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors