

### Donkey Anti-Sheep IgG H&L (Texas Red ®) ab6898

#### 概述

产品名称	驴抗羊 IgG H&L (Texas Red ®)
宿主	Donkey
靶标种属	Sheep
经测试应用	适用于: ICC/IF, IHC-P, IHC-Fr, Immunomicroscopy, Flow Cyt, ELISA
免疫原	Full length native Sheep IgG (purified).
偶联物	Texas Red ®. Ex: 596nm, Em: 620nm

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 6.50 Preservative: 0.01% Sodium azide Constituents: 0.42% Tripotassium orthophosphate, 0.87% Sodium chloride, 1% BSA
纯度	Affinity purified
纯化说明	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Sheep IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities.
共轭说明	Texas Red™ Sulfonyl Chloride (Molecular Weight 625 daltons) Absorption Wavelength: 596 nm Emission Wavelength: 620 nm Fluorochrome/Protein Ratio: 4.9 moles Texas Redä per mole of Rabbit IgG
克隆	多克隆

#### 应用

##### The Abpromise guarantee

**Abpromise™** 承诺保证使用 ab6898 于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		
IHC-P		
IHC-Fr		
Immunomicroscopy		
Flow Cyt		
ELISA		

#### 应用说明

Flow Cyt: Use at an assay dependant concentration.

IHC-P: Use at an assay dependant concentration.

IHC-Fr: Use at an assay dependant concentration.

IM: Use at an assay dependant concentration.

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

#### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors