

Goat Anti-Armenian hamster IgG H&L (Biotin) ab5744

9 References 1 图像


概述	
产品名称	山羊抗亚美尼亚仓鼠IgG H&L (Biotin)
宿主	Goat
靶标种属	Armenian hamster
特异性	Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-biotin, anti-Goat Serum, Armenian Hamster IgG and Armenian Hamster Serum. Cross reactivity against mouse has not been determined.
经测试应用	适用于: WB, ICC/IF, Dot blot, ELISA, Immunomicroscopy, IHC-P
免疫原	Armenian Hamster IgG whole molecule
偶联物	Biotin
性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 6.50 Preservative: 0.01% Sodium azide Constituents: 0.42% Tripotassium orthophosphate, 0.87% Sodium chloride, 1% BSA
纯化说明	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Armenian Hamster IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities.
共轭说明	Biotinamidocaproate N-Hydroxysuccinimide Ester (BAC) Biotin/Protein Ratio: 10-20 BAC molecules per Goat IgG molecule.
克隆	多克隆
同种型	IgG
应用	

The Abpromise guarantee **Abpromise™**承诺保证使用ab5744于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/2000 - 1/10000.
ICC/IF		1/1000 - 1/5000.
Dot blot		Use at an assay dependent dilution.
ELISA		1/250000.
Immunomicroscopy		Use at an assay dependent dilution.
IHC-P		Use at an assay dependent dilution.

图片



ab5744 was used at dilution 1/500 with the primary antibody **ab33483** in IHC-P. See the review on **ab33483**.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Goat Anti-Armenian hamster IgG H&L (Biotin) (ab5744)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors