

Goat Anti-Chicken IgY H&L (FITC) ab46969

[6 References](#) [1 图像](#)

概述

产品名称	山羊抗鸡IgY H&L (FITC)
宿主	Goat
靶标种属	Chicken
经测试应用	适用于: ICC/IF, IHC-Fr
免疫原	Chicken IgY
偶联物	FITC. Ex: 493nm, Em: 528nm

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7 Preservative: 0.02% Sodium azide Constituents: PBS, 0.2% BSA
纯度	Isotonic 0.9%, w/v)
纯化说明	Affinity purified
克隆	Ab46969 is specific to chicken IgY and will not cross react with mammalian IgG.
同种型	多克隆 IgG

应用

The Abpromise guarantee **Abpromise™** 承诺保证使用ab46969于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		

应用	Ab评论	说明
IHC-Fr		

应用说明


ICC/IF: 1/250.

IHC-Fr: Use at an assay dependent dilution (PMID 18583367).

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

图片



ab4674 staining GFAP in human brain cells by Flow Cytometry.

Cells were cultured in stem cell media RHB-A and collected using 0.05% Trypsin. The sample was incubated with the primary antibody (1/100 in PBS + 10% FCS) for 20 minutes at 25°C. A FITC-conjugated Goat anti-chicken IgY (ab46969) (1/100) was used as the secondary antibody.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors