

Goat Anti-Chicken IgY H&L (FITC) ab6873

★★★★★ [1 Abreviews](#) [5 References](#) [2 图像](#)

概述

产品名称	山羊抗鸡IgY H&L (FITC)
宿主	Goat
靶标种属	Chicken
经测试应用	适用于: Immunomicroscopy, Flow Cyt, ICC/IF, IHC-P, IHC-Fr, ELISA
免疫原	chicken IgY whole molecule
偶联物	FITC. Ex: 493nm, Em: 528nm

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 6.50 Preservative: 0.01% Sodium azide Constituents: 0.88% Sodium chloride, 1% BSA, 0.42% Tripotassium orthophosphate
纯度	Affinity purified
纯化说明	This product was prepared from monospecific antiserum by immunoaffinity chromatography using chicken IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities.
共轭说明	Diaminotriazinylaminofluorescein (DTAF) (Molecular Weight 530 daltons) Absorption Wavelength: 495 nm Emission Wavelength: 528 nm Fluorochrome/Protein Ratio: 2.3 moles DTAF per mole of Goat IgG
克隆	多克隆


应用

The Abpromise guarantee **Abpromise™** 承诺保证使用ab6873于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
Immunomicroscopy		Use at an assay dependent concentration.
Flow Cyt		Use at an assay dependent dilution.
ICC/IF	★★★★★ (1)	1/1000 - 1/5000.
IHC-P		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.
ELISA		1/10000 - 1/50000.

图片



Immunocytochemistry/ Immunofluorescence - Goat
Anti-Chicken IgY H&L (FITC) (ab6873)

Ab6873 was used at dilution 1/500 with the primary antibody [ab39376](#) in ICC/IF. See the review on [ab39376](#).



Immunohistochemistry (Frozen sections) - Goat
Anti-Chicken IgY H&L (FITC) (ab6873)

Ab6873 was used at dilution 1/500 with the primary antibody [ab13991](#) in IHC-Fr. See the review on [ab13991](#).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors