

Goat Anti-Cow IgG H&L ab102150

1 References

| | |
|-------|---|
| 概述 | |
| 产品名称 | 山羊抗牛IgG H&L |
| 宿主 | Goat |
| 靶标种属 | Cow |
| 特异性 | By immunoelectrophoresis and ELISA this antibody reacts specifically with bovine IgG and with light chains common to other bovine immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. This antibody may cross react with IgG from other species. |
| 经测试应用 | 适用于: ICC, IHC-P, ELISA, WB |

| | |
|------|--|
| 性能 | |
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. |
| 存储溶液 | pH: 6.8 Preservative: 0.09% Sodium azide Constituents: 0.2% BSA, PBS |
| 纯度 | Immunogen affinity purified |
| 纯化说明 | This antibody was isolated by affinity chromatography using antigen coupled to agarose beads. |
| 克隆 | 多克隆 |
| 同种型 | IgG |

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab102150于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

| 应用 | Ab评论 | 说明 |
|-------|------|-----------------|
| ICC | | 1/200 - 1/2000. |
| IHC-P | | 1/200 - 1/2000. |

| 应用 | Ab评论 | 说明 |
|-------|------|---|
| ELISA | | 1/100 - 1/30000. (Coating: 1/100 - 1/500; Primary: 1/1000 - 1/30000). |
| WB | | 1/1000 - 1/30000. (Colorimetric: 1/1000 - 1/10000; Chemiluminescent: 1/1000 - 1/30000). |

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors