

Goat Anti-Cow IgG H&L preadsorbed ab102159

概述

产品名称	山羊抗牛IgG H&L预吸附二抗
宿主	Goat
靶标种属	Cow
特异性	By immunoelectrophoresis and ELISA this antibody reacts specifically with bovine IgG and with light chains common to other bovine immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. Reduced cross-reactivity to human, mouse and rat IgG was detected. This antibody may cross react with IgG from other species.
经测试应用	适用于: ELISA, IHC-P, ICC, WB
最小交叉反应	Human, Mouse, Rat

[more details](#)

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 6.8 Preservative: 0.09% Sodium azide Constituent: PBS
纯度	Immunogen affinity purified
纯化说明	Antiserum was cross adsorbed using human, mouse and rat immunosorbents to remove cross reactive antibodies. This antibody was isolated by affinity chromatography using antigen coupled to agarose beads.
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab102159于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		1/100 - 1/30000. Coating 1/100 - 1/500. Primary 1/1000 - 1/30000
IHC-P		1/200 - 1/2000.
ICC		1/200 - 1/2000.
WB		1/1000 - 1/30000. Colorimetric 1/1000 - 1/10000. Chemiluminescence 1/1000 - 1/30000.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors