

Goat F(ab')<sub>2</sub> Anti-Rabbit IgG - H&L (AP), pre-adsorbed ab98505

6 References

概述	
产品名称	山羊F(ab') <sub>2</sub> 抗兔IgG - H&L (AP),预吸附二抗
宿主	Goat
靶标种属	Rabbit
经测试应用	适用于: ICC, IHC-P, ELISA, WB
最小交叉反应	Chicken, Cow, Goat, Horse, Human, Mouse, Rat, Sheep
偶联物	Alkaline Phosphatase
<a href="#">more details</a>	

性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: 0.2% BSA, 0.58% Sodium chloride, 0.0095% Magnesium chloride, 0.00136% Zinc chloride, 1.19% HEPES
纯度	Immunogen affinity purified
纯化说明	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads. F(ab') <sub>2</sub> fragment were generated using a pepsin digestion. Fc fragments and whole IgG molecules have been removed. Fragments were conjugated to Alkaline Phosphatase.
克隆	多克隆
同种型	IgG

应用	
The Abpromise guarantee	<u>Abpromise™</u> 承诺保证使用ab98505于以下的经测试应用
“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。	

应用	Ab评论	说明
ICC		Use at an assay dependent dilution.
IHC-P		1/100 - 1/1000.
ELISA		1/1000 - 1/10000. (Primary ELISA)
WB		1/1000 - 1/10000. (Chemiluminescent)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors