

Goat Anti-Ferret IgG H&L (FITC) ab112769

2 References

概述	
产品名称	山羊Anti-Ferret IgG H&L (FITC)
宿主	Goat
靶标种属	Ferret
特异性	By immunoelectrophoresis and ELISA ab112769 reacts specifically with Ferret IgG and with light chains common to other Ferret immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. May cross react with IgG from other species.
经测试应用	适用于: IHC-P, ICC/IF
免疫原	Other Immunogen Type corresponding to IgG.
偶联物	FITC. Ex: 493nm, Em: 528nm
性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: 0.2% BSA, 99% PBS
纯度	Immunogen affinity purified
纯化说明	ab112769 was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to fluorescein isothiocyanate (FITC).
共轭说明	F/P ratio is 3 to 7.
克隆	多克隆
同种型	IgG
应用	

The Abpromise guarantee **Abpromise™** 承诺保证使用ab112769于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P		1/50 - 1/500.
ICC/IF		1/50 - 1/500.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors