



Product datasheet

Goat Anti-Human IgA (SureLight® PE) ab99909

1 References

概述

产品名称	山羊抗人IgA (SureLight® PE)
宿主	Goat
靶标种属	Human
经测试应用	适用于: Flow Cyt
偶联物	SureLight® PE

性能

形式	Lyophilized:Lyophilized material is stable for up to one year. Reconstitue to 1.0 ml with distilled deionized water vortex gently and allow vial to sit on ice for 20 minutes. After product has been reconstituted, product should be stored at 2-8°C in the dark and be used within 3 month. If further dilution of the conjugate is required, use diluted material within one week.
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.20 Preservative: 0.05% Sodium azide Constituents: 0.1% BSA, 0.58% Sodium chloride, 1.64% Sodium phosphate, 1.71% Sucrose
纯度	Immunogen affinity purified
共轭说明	SureLight® R-Phycoerythrin Excitation max: 565>498 nm Emission max.: 578 nm F/P 1.3:1
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab99909于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		1/100 - 1/1000.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors