

Goat Anti-Human IgE epsilon chain ab9159

5 References

概述

产品名称	山羊抗人IgE epsilon chain
宿主	Goat
靶标种属	Human
特异性	The affinity purified anti IgE reacts only with IgE (several different IgE kappa and lambda myeloma proteins) by IEP and gel diffusion techniques.
经测试应用	适用于: ELISA, Immunodiffusion, Immunoelectrophoresis, WB
免疫原	Goats were immunized with a mixture of several highly purified human IgE myeloma proteins in Freund's adjuvant

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.09% Sodium azide
纯度	Immunogen affinity purified
纯化说明	Serum was collected and passed through several non-IgE containing immunoaffinity supports to make the preparation delta chain specific. The anti delta was immunoaffinity purified off an immunoabsorbant containing a different IgE than was used for immunization thus minimizing idiotypic influences.
克隆	多克隆

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab9159于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		Use at an assay dependent concentration.

应用	Ab评论	说明
Immunodiffusion		Use at an assay dependent concentration.
Immunoelectrophoresis		Use at an assay dependent concentration.
WB		Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors