

# Goat Anti-Human IgE H&L (HRP) ab73901

## 3 References

### 概述

产品名称	山羊抗人IgE H&L (HRP)
宿主	Goat
靶标种属	Human
经测试应用	适用于: ELISA
免疫原	Full length protein, highly purified human IgE.
偶联物	HRP

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Constituents: 0.2% BSA, 99.75% PBS, 0.05% CMIT/MIT based preservative
纯度	Protein G purified
克隆	多克隆
同种型	IgG

### 应用

**The Abpromise guarantee**      **Abpromise™**承诺保证使用ab73901于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		1/10000 - 1/100000.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors