

Goat Anti-Human IgE (SureLight® APC) ab99898

1 图像

概述

产品名称	山羊抗人IgE (SureLight® APC)
宿主	Goat
靶标种属	Human
经测试应用	适用于: Flow Cyt
免疫原	Full length native protein (purified) corresponding to IgE.
偶联物	SureLight® APC. Ex: 652nm, Em: 657nm

性能

形式	Lyophilized:Reconstitute with 1 ml of ddH2O
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 8.20 Constituents: 1.71% Sucrose, 0.1% BSA, 0.87% Sodium chloride, 0.121% Tris
纯度	Immunogen affinity purified
共轭说明	Fluor:Protein ~2:1 Visible absorption maxima 652 Emission maximum 657
克隆	多克隆
同种型	IgG
常规说明	After product has been reconstituted, product should be stored at 2-8°C in the dark and be used within 3 months. If further dilution of the conjugate is required, use diluted material within one week.


应用

The Abpromise guarantee **Abpromise™** 承诺保证使用ab99898于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use a concentration of 3 µg/ml.

图片



Flow Cytometry - Goat Anti-Human IgE (SureLight® APC) (ab99898)

Fluorescence excitation and emission spectra of Allophycocyanin in 100mM sodium phosphate (pH7.2) + 1mM EDTA and 1mM sodium azide. The emission scan was taken with excitation at 630nm. The excitation scan was taken with an emission at 660nm. The curves were normalized to equalize peak heights.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors