

Goat Anti-Human IgG (Fab')₂ (HRP) ab87422

[6 References](#) [1 图像](#)

概述

产品名称	山羊抗人IgG (Fab') ₂ (HRP)
宿主	Goat
靶标种属	Human
特异性	ab87422 is specific for the Fab and F(ab') ₂ portion of human IgG by IEP. It reacts with other human immunoglobulins through common light chain reactivity.
经测试应用	适用于: ELISA, WB
免疫原	Other Immunogen Type corresponding to IgG. Human serum.
偶联物	HRP

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.01% Thimerosal (merthiolate) Constituents: 0.1% BSA, PBS
纯度	Protein G purified
纯化说明	This antibody has been cross adsorbed against the Fc domain and will react with less than 1 % of the Fc domain of human IgG.
克隆	多克隆
同种型	IgG

应用


The Abpromise guarantee **Abpromise™** 承诺保证使用ab87422于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		1/1000 - 1/10000.

应用	Ab评论	说明
WB		1/1000 - 1/10000.

图片



Characterisation of virus activity and functional expression of monoclonal antibodies.

(Panel D) IgG1 antibody expression, assessed by anti-IgG1 ELISA, in cellular supernatants of carcinoma cell lines infected for 24-96hrs with NG-135 at 1ppc.

Media was removed from infected or control cells, clarified by centrifugation and diluted 1 in 2 in 3% BSA/PBS for storage at -20°C. For quantification of human IgG1, 96 well plates were coated overnight, 4°C, with mouse monoclonal anti-human IgG1 Fc antibody (**ab1927**, Abcam) at a 1:1000 dilution in carbonate/bicarbonate buffer. Plates were blocked in 3% BSA/PBS, washed and samples and standards added to the plate for 1 hr and RT. Secondary detection was carried out using a HRP conjugated goat anti-human IgG Fab (ab87422, Abcam) incubated for 1 hr, RT.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors