

Goat Anti-Human IgG Fc (DyLight® 650) preadsorbed ab98622

8 References 1 图像

概述		
产品名称	山羊抗人IgG Fc (DyLight® 650)预吸附二抗	
宿主	Goat	
靶标种属	Human	
经测试应用	适用于: IHC-P, ICC/IF, Flow Cyt	
最小交叉反应	Chicken, Cow, Goat, Horse, Mouse, Rabbit, Rat, Sheep	
偶联物	DyLight® 650. Ex: 654nm, Em: 673nm	
		more details

性能		
形式	Liquid	
存放说明	Shipped at 4°C. Store at +4°C.	
存储溶液	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: PBS, 0.2% BSA	
纯度	Immunogen affinity purified	
纯化说明	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to DyLight® 650.	
克隆	多克隆	
同种型	IgG	
常规说明	DyLight® 650 replaces DyLight® 649.	

应用


The Abpromise guarantee **Abpromise™**承诺保证使用ab98622于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明

应用	Ab评论	说明
IHC-P		1/50 - 1/500. DyLight® 650 replaces DyLight® 649. The DyLight® 650 excitation is 652nm and the emission is 672nm. DyLight® 650 hence provides the same far-red fluorescence and photostability as the DyLight® 649 dye. Please refer to the vial in order to check whether the product you received is DyLight® 649 or DyLight® 650. This information however will have no impact on your
ICC/IF		1/50 - 1/500.
Flow Cyt		1/50 - 1/200.

图片



Emission spectra of DyLight® fluorochromes available in our catalog.

Line colors represent the approximate visible colors of the wavelength maxima.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors